

# Rotaid WALL CABINETS FOR THE HeartSave FAMILY


## AED WALL CABINETS Rotaid PLUS FOR INDOOR USE



The Rotaid can be opened easily, quickly and intuitively with a simple rotary motion. The cover is tethered to the back of the Rotaid and therefore will not get in the way in case the device has to be used. The first responder immediately has complete access to the AED.

### FEATURES Rotaid PLUS

- ▶ 100% polycarbonate with a thickness of 5 mm for solid protection from water and dust
- ▶ A cover with UV filter provides protection from direct and indirect sunlight.
- ▶ The Rotaid Plus can optionally be equipped with a visible and acoustic alarm (85 dB)
- ▶ Dimensions: 47 cm x 47 cm x 18 cm
- ▶ Weight: 2.1 kg
- ▶ Protection class IK10 (Impact)



# AED WALL CABINET Rotaid SOLID PLUS FOR OUTDOOR USE

## FEATURES Rotaid SOLID PLUS

- ▶ Equipped with a heating device: The heating device ensures that the temperature on the inside of the Rotaid is at least 5°C, even if the outside temperature drops to -30° .
- ▶ 100% polycarbonate with a thickness of 5 mm for solid protection from water and dust
- ▶ A cover with UV filter provides protection from direct and indirect sunlight.
- ▶ The Rotaid SOLID PLUS is equipped with a visible and acoustic alarm (85 dB)
- ▶ Dimensions: 47cm x 47cm x 18cm
- ▶ Weight: 7.1 kg
- ▶ Protection class: IP54
- ▶ Protection class IK10 (Impact)



## Rotaid PLUS OPEN



**METRAX GmbH**  
Rheinwaldstr. 22  
D-78628 Rottweil  
Germany

Phone: +49 741 257-0  
Fax: +49 741 257-235  
info@primedic.com

[www.primedic.com](http://www.primedic.com)



**PRIMEDIC**  
Saves Life. Everywhere.